



# Recombinant Human UPF0464 protein C15orf44 (C15orf44)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-YP847227HU  |
| <b>Abbreviation</b>      | C15orf44  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q96SY0  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MPTVVVMDVS LSMTRPVSIE GSEEYQRKHL AAHGLTMLFE HMATNYKLEF<br>TALVVFSSLW ELMVPFTRDY NTLQEALSNM DDYDKTCLES ALVGVCNIVQ<br>QEWGGAIPCQ VVLVTDGCLG IGRGSLRHSL ATQNQRSESN RFPLPFPFSP<br>KLYIMCMANL EELQSTDSLE CLERLIDLNN GEGQIFTIDG PLCLKNVQSM<br>FGKLIDLAYT PFHAVLKCGRH LTADVQVFPR PEPFVVDEEI DPIPKVINTD<br>LEIVGFIDIA DISSPPVLSR HLVLPIALNK EGDEVGTGIT DDNEDENSAN<br>QIAGKIPNFC VLLHGSLKVE GMVAIVQLGP EWHGMLYSQA DSKKKSNLMM<br>SLFEPGPEPL PWLGKMAQLG PISDAKENPY GEDDNKSPFP LQPKNKRSYA<br>QNVTVWIKPS GLQTDVQKIL RNARKLPEKT QTFYKELNRL RKAALAFGFL<br>DLLKGVADML ERECTLLPET AHPDAAFQLT HAAQQLKLAS TGTSEYAAAYD<br>QNITPLHTDF SGSSTERI |
| <b>Source</b>            | Yeast   |
| <b>Target Names</b>      | INTS14  |
| <b>Protein Names</b>     | Recommended name: UPF0464 protein C15orf44  |
| <b>Expression Region</b> | 1-518   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life   |



of lyophilized form is 12 months at -20°C/-80°C.