



Recombinant Hypocrea jecorina Transcription factor ACEII (ace2)

Product Code	CSB-MP847239HYE
Abbreviation	ace2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96WN6
Product Type	Recombinant Protein
Immunogen Species	Hypocrea jecorina (Trichoderma reesei)
Purity	>85% (SDS-PAGE)
Sequence	MDLRQACDRC HDKKLRCPRI SGSPCCSRCA KANVACVFSP PSRPFRPHEP LNHSHEHSHS HSHNHNGVGV SFDWLDLMSL EQQQEQQQGQ PQHPPPPVQT LSERLAALLC ALDRMLQAVP SSLDMHHVSR QQLREYADTV GTGFDLQSTL DSLLHHAQDL ASLYSEAVPA SFNKRTTAAE ADALCAVPDC VHQDRTSLHT TPLPKLDHAL LNLVMACHIR LLDVMDTLAE HGRMCAFMVA TLPPDYDPKF AVPEIRVGTF VAPTDTAASM LLSVVVELQT VLVARVKDLV AMVDQVKDDA RAAREAKVVR LQCGILLERA ESTLGEWSRF KDGLVSARLL K
Source	Mammalian cell
Target Names	ace2
Protein Names	Recommended name: Transcription factor ACEII
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.