



Recombinant Arabidopsis thaliana Actin-depolymerizing factor 12 (ADF12)

Product Code	CSB-MP847516DOA
Abbreviation	ADF12
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8LFH6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MANAASGMAV EDECKLKFLE LKAKRNYRFI IFRIDGQQVV VEKLGSPQEN YDDFTNYLPP NECRYAVYDF DFTTAENIQK SKIFFIAWSP DSSRVRMKMV YASSKDRFKR ELDGIQVELQ ATDPSEMSLD IIKSRAL
Source	Mammalian cell
Target Names	ADF12
Protein Names	Recommended name: Actin-depolymerizing factor 12 Short name= ADF-12 Short name= AtADF12
Expression Region	1-137
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.