



# Recombinant Human Uncharacterized protein C1orf87 (C1orf87)

<b>Product Code</b>	CSB-EP847603HU
<b>Abbreviation</b>	C1orf87
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N0U7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSSAWKTPRG SDAMPEIMVK IIGSKHFQYL VEKPKIKEND SLKTETQTMH QKPMTDNARQ MSRDTPVPIN FTDQQTDDNP DDVKEKKHPE NNQKSENNQK LLTGANSSRF LDGNIPSQAN VHCSSVPTGD QSLSYVHGIP RRKLRDWSLE QMVRGSSDQP EDIGQSPSGT TNEDAFLLAL VRRELKSRPL SSNLLEKLQK ELKILDPISS GFLLQSQLSR LFLKHEVPLQ LPTVKILCQR FSKRGSPEMV NYEKLLWFLN SAASDYPQQN KAAADLRKTE SHGTHSQSTP PQHSSSQPEV NRSLLLEILKM ALRTTNGRLN IDNLNLSFRK EDRSFSGCLP LPKVRAICGK HGLYLTLSELL ETLLNHQDLG YQNEIKWQNF VEMLTRASSD LLSDLPTGKN EKKAPAPPME PEVPEMSQSK TEHMKTPEEE LQPESSPAET SACKDPLKPL KIRPVSQPFV NPAVKNKAAE CETWIDRFRK LENALYLCDL SNTGVLEKER ARRLIHNYNL IYNLSLSPQK IDQALRRFRS GENMLLEPAL RYLKEL
<b>Source</b>	E.coli
<b>Target Names</b>	C1orf87
<b>Protein Names</b>	Recommended name: Uncharacterized protein C1orf87
<b>Expression Region</b>	1-546
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.