



Recombinant Human Beckwith-Wiedemann syndrome chromosomal region 1 candidate gene B protein (SLC22A18AS)

Product Code	CSB-EP847610HU-B
Abbreviation	SLC22A18AS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N1D0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGELPGSEGM WENCPLGWVK KKASGTLAPL DFLLQRKRLW LWASEPVRPQ PQGIHRFREA RRQFCRMRGS RLTGGRKGF SSGLRFGRGG FSEEVMPQPV LKAMRCAEGA WWFSPDGPAG SAASIWPAEG AEGLPGQLGR DRLEVVYSVP DNVPQNGSR RPLVCKITGK CL SVCSEENA KAGGCSAFPL LLSQLGARMT GREHAHKGPE LTT PDSGLPR PPNPALAGFR ALAQHSPPLG TSTPSAVLLS AAT
Source	E.coli
Target Names	SLC22A18AS
Protein Names	Recommended name: Beckwith-Wiedemann syndrome chromosomal region 1 candidate gene B protein Alternative name(s): Organic cation transporter-like protein 2 antisense protein Solute carrier family 22 member 1-like antisense protein Solute ca
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.