



Recombinant Human E3 ubiquitin-protein ligase pellino homolog 3 (PELI3)

Product Code	CSB-EP847618HU-B
Abbreviation	PELI3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N2H9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVLEGNPEVG SPRTSDLQHR GNKGSCVLSS PGEDAQPGEE PIKYGELIVL GCCEEGGEET EAQRGEVTGP RAHSCYNGCL ASGDKGRRRS RLALSRRSHA NGVKPDVMHH ISTPLVSKAL SNRGQHSISY TLSRSHSVIV EYTHDSDTDM FQIGRSTENM IDVVVTDTSPP GGGAAEGPSA QSTISRYACR ILCDRRPPYT ARIYAAGFDA SSNIFLGERA AKWRTPDGLM DGLTTNGVLV MHPAGGFSED SAPGVWREIS VCGNVYTLRD SRSAQQRGKL VENESNVLQD GSLIDLGCAT LLWRTPAGLL RAPTLKQLEA QRQEANAARP QCPVGLSTLA FPSPARGRTA PDKQQPWVYV RCGHVHGYHG WGCRRERGPQ ERECPLCRLV GPYVPLWLGQ EAGLCLDPGP PSHAFAPCGH VCSEKTARYW AQTPLPHGTH AFHAACPFCC AWLTGEHGCV RLIFQGPLD
Source	E.coli
Target Names	PELI3
Protein Names	Recommended name: E3 ubiquitin-protein ligase pellino homolog 3 Short name= Pellino-3 EC= 6.3.2.-
Expression Region	1-469
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.