



Recombinant Human MAGUK p55 subfamily member 5 (MPP5)

Product Code	CSB-BP847626HU
Abbreviation	MPP5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N3R9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MTTSHMNGHV TEESDSEVKN VDLASPEEHQ KHREMAVDCP GDLGTRMMPI RRSAQLERIR QQQEDMRRRR EEEGKKQELD LNSSMRLKKL AQIPPKTGID NPMFDTEEGI VLESPHYAVK ILEIEDLFSS LKHIQHTLVD SQSQEDISLL LQLVQNKDFQ NAFKIHNAIT VHMNKASPPF PLISNAQDLA QEVQTVLKPV HHKEGQELTA LLNTPHIQAL LLAHDKVAEQ EMQLEPITDE RYYESIGQYG GETVKIVRIE KARDIPLGAT VRNEMDSVII SRIVKGGAAE KSGLLHEGDE VLEINGIEIR GKDVNEVFDL LSDMHGTLTF VLIPSQQIKP PPAKETVIHV KAHFDYDPSD DPYVPCRELG LSFQKGDILH VISQEDPNWW QAYREGDEDN QPLAGLVPGK SFQQQREAMK QTIEEDKEPE KSGKLWCAKK NKKKRKKVLY NANKNDDYDN EEILTYEEMS LYHQPANRKR PIILIGPQNC GQNELRQRLM NKEKDRFASA VPHTTRSRRD QEVAGRDIYHF VSRQAFEADI AAGKFIEHGE FEKNLYGTSI DSVRQVINSK KICLLSLRTQ SLKTLRNSDL KPYIIFIAPP SQERLRALLA KEGKNPKPEE LREIIEKTRE MEQNNGHYFD TAIVNSDLDK AYQELLRLIN KLDTEPQWVP STWLR</p>
Source	Baculovirus
Target Names	MPP5
Protein Names	Recommended name: MAGUK p55 subfamily member 5
Expression Region	1-675
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.