



Recombinant Human F-box/WD repeat-containing protein 8 (FBXW8)

Product Code	CSB-MP847627HU
Abbreviation	FBXW8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N3Y1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDDYSLDEFRRRWQEELAQA QAPKKRRRPE AAERRARRPE VGSGRGEQAS GDPALAQRLLEGAGRPPAAR ATRAEGQDVA SRSRSPLARE GAGGGEQLVD QLIRDLNEMN DVPFFDIQLP YELAINIFQY LDRKELGRCA QVSKTWK VIA EDEVLYWRLC QQEGHLPDSS ISDYSCWKLI FQECRAKEHM LRTNWKNRKG AVSELEHVPD TVLCDVHSHD GVVIAGYTSG DVRVWDTRTW DYVAPFLESE DEEDEPGMQP NVSFVRINSS LAVAAYEDGF LNIWDLRTGK YPVHRFEHDA RIQALALSQD DATVATASAF DVVMLSPNEE GYWQIAAEFE VPKLVQYLEI VPETRRYPVA VAAAGDLMYL LKAEDSARTL LYAHGPPVTC LDVSANQVAF GVQGLGWVYE GSKILVYSLE AGRRLKLG VLRDFTCVNL SDSPNLMVS GNMDGRVRIH DLRSGNIALS LSAHQLRVSA VQMDDWKIVS GGEEGLVSVW DYRMNQKLWE VYSGHPVQHI SFSSHSLITA NVPYQTVMRN ADLDSFTTHR RHRGLIRAYE FAVDQLAFQS PLPVCRSSCD AMATHYYDLA LAFPYNHV
Source	Mammalian cell
Target Names	FBXW8
Protein Names	Recommended name: F-box/WD repeat-containing protein 8 Alternative name(s): F-box and WD-40 domain-containing protein 8 F-box only protein 29
Expression Region	1-598
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the F-box protein family, members of which are characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into three classes: Fbws



containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. This protein contains a WD-40 domain, in addition to an F-box motif, so it belongs to the Fbw class. Alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.