



Recombinant Human Uncharacterized protein C4orf37 (C4orf37)

Product Code	CSB-BP847628HU
Abbreviation	C4orf37
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N412
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MYDRAPRLK LAEGGSTEAH VGPGSYQVPF LKQQATGSNA PFLSLTARES TFTIASSIEK AVPGPGHYNV SEAQKISRSP TLTRSVDVPS IPSCGKSYGY HINDDGSIIK CFPPACDSTL GPAYYKQFD VSNATLKYKG IHFGNSSGRQ ELPKKSGPGP GQYDIVQKKT SYYENVNIKR DQQQNYCSFI PRLYEIIVLQ EKKKRFLPMK SITPAPGTYN EPRTALKSLK KTSGLKNIPF GQSAVRFTQD IRTEEMPGPG FYNVLNNTII ASVRNICKSK QKKSFAFGSSV PRTFFSVQKE ACATPGPADY QEFWHSQGVG ISDELPNLTN KYAAFLSRAK RTMKVPDMVI PAPGSYDVHK SYEMSQVKHK YMPRSLVAK RKHASFLSAT PRCLEKVTGD PGPAAYNPVL RKSCPIPLFV KASKRFEEK EITPGPATYE ISQEKKKGNL IGEMAADIM
Source	Baculovirus
Target Names	STPG2
Protein Names	Recommended name: Uncharacterized protein C4orf37
Expression Region	1-459
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.