



Recombinant Human Spermatogenesis associated 6-like protein (SPATA6L)

Product Code	CSB-EP847631HU-B
Abbreviation	SPATA6L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N4H0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPLEVVVELQ IRAISCPGVF LPGKQDVYLG VYLMNQYLET NSFPSAFPIM IQESMRFEKV FESAVIDPGAV VDLLEMWDEL AYYEENTRDF LFPEPKLTPS HPRRCREVLK KTAGLFGPIA PKIEFSTRTA IRECVFLHRN RFLEERHESR RPLSTSHPI FPLNTIKMKL KENNLNRLPK GMQARAPSQY STRHFFQDQP AQLNLGNFVK ISGGSKPPFV VRHVDSAKPF GENISEHHLR RSRKSKFSD FPFPTRRASS LDSLAANVKV IKEPDERIVL RSDSSSCLDS SQFGKSSSSK QGDADFHGKA SFATYQHSTS PGPLDQPLL RERHFGSQST WKNIHERVCS LLTSHRAQLH QNKEDSTSEV NYIIRPSYP LKKYSLHEQR YF
Source	E.coli
Target Names	SPATA6L
Protein Names	Recommended name: Spermatogenesis associated 6-like protein
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.