



# Recombinant Human Coiled-coil domain-containing protein 82 (CCDC82)

<b>Product Code</b>	CSB-YP847633HU
<b>Abbreviation</b>	CCDC82
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N4S0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MIHVRRHETR RNSKSHVPEQ KSRVDWRRTK RSSISQLLDS DEELDSEEFD SDEELSDDES FENDEELDSN KGPDCNKTPG SERELNLSKI QSEGNSKCL INSGNGSTYE EETNKIKHRN IDLQDQEKHL SQEDNDLNKQ TGQIIEDDQE KHLSEQEDNDL NKQTGQIIED DLEEEDIKRG KRKRLSSVMC DSDESDDSDI LVRKVGVKRP RRVVEDEGSS VEMEQKTPEK TLAQKREKL QKLKELSKQR SRQRRSSGRD FEDSEKESCP SSDEVDEEEE EDNYEDEDG DDYIIDDFV QDEEGDEENK NQQGEKLTTS QLKLVKQNSL YSFSDHYTHF ERVVKALLIN ALDESFLGTL YDGTRQKSYA KDMLTSLHYL DNRVQPRLE SLVSRSRWKE QYKERVERNYS NVSIHLKNPE NCSCQACGLH RYCKYSVHLS GELYNTRTMQ IDNFMSHDKQ VFTVGRICAS RTRIIYHKLKH FKFELYQECC TIAMTEEVED EQVKETVERI FRRSKENGWI KEKYGQLEEY LNFADYFQEE KFEI
<b>Source</b>	Yeast
<b>Target Names</b>	CCDC82
<b>Protein Names</b>	Recommended name: Coiled-coil domain-containing protein 82
<b>Expression Region</b>	1-544
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.