



Recombinant Human T-complex protein 11-like protein 2 (TCP11L2)

Product Code	CSB-MP847634HU
Abbreviation	TCP11L2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N4U5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPFNGEKQCV GEDQPSDSDS SRFSESMASL SDYECSRQSF ASDSSSKSSS PASTSPPRVV TFDEVMATAR NLSNLTLAHE IAVNENFQLK QEALPEKSLA GRVKHIVHQA FWDVLDSELN ADPPEFEHAI KLFEEIREIL LSFLTPGGNR LRNQICEVLD TDLIRQQA EH SAVDIQGLAN YVISTMGKLC APVRDNDIRE LKATGNIVEV LRQIFHVLDL MQMDMANFTI MSLRPHLQRR LVEYERTKFQ EILEETPSAL DQTTEWIKES VNEELFSLSE SALTPGAENT SKPSLSPTLV LNNSYLKLLQ WDYQKKELPE TLMTDGARLQ ELTEKLNQLK IIACLSLITN NMVGAITGGL PELASRLTRI SAVLLEGMNK ETFNLKEVLN SIGIQCVEV NKTLMERGLP TLNAEIQANL IGQFSSIEEE DNPIWSLIDK RIKLYMRLL CLPSPQKCM PMPGGLAVIQ QEALGSQY ANIVNLNKQV YGPFIYANILR KLLFNEEAMG KVDASPTN
Source	Mammalian cell
Target Names	TCP11L2
Protein Names	Recommended name: T-complex protein 11-like protein 2
Expression Region	1-519
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.