



# Recombinant Human 2-oxoglutarate and iron-dependent oxygenase domain-containing protein 1 (OGFOD1)

<b>Product Code</b>	CSB-BP847635HU
<b>Abbreviation</b>	OGFOD1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N543
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNGKRPAEPG PARVGKKGKK EVMAEFSDAV TEETLKKQVA EAWSRRTPFSS HEVIVMDMDP FLHCVIPNFI QSQDFLEGLQ KELMNLDFHE KYNDLYKFQQ SDDLKKRREP HISTLRKILF EDFRSWLSDI SKIDLESTID MSCAKYEFTD ALLCHDDELE GRRIAFILYL VPPWDRSMGG TLDLYSIDEH FQPKQIVKSL IPSWNKL VFF EVSPVSFHQV SEVLSEEKSR LSISGWFHGP SLTRPPNYFE PPIPRSPHIP QDHEILYDWI NPTYLDMDYQ VQIQEEFEES SEILLKEFLK PEKFTKVCEA LEHGHVEWSS RGPPNKR FYE KAEESKLPEI LKECMKLFRS EALFLLLSNF TGLKLHFLAP SEEDEMNDKK EAETTDITEE GTSHSPPEPE NNQMAISNNS QQSNEQTDPE PEENETKKES SVPMCQGELR HWKTGHYTLI HDHСКАEFAL DLILYCGCEG WEPEYGGFTS YIAKGEDEEL LTVNPESNSL ALVYRDRETL KFKVHINHRS LEQKKTFFNR TGFWDFFSIY YE
<b>Source</b>	Baculovirus
<b>Target Names</b>	OGFOD1
<b>Protein Names</b>	Recommended name: 2-oxoglutarate and iron-dependent oxygenase domain-containing protein 1 EC= 1.14.11.- Alternative name(s): Termination and polyadenylation 1 homolog
<b>Expression Region</b>	1-542
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.