



Recombinant Human 2-oxoglutarate and iron-dependent oxygenase domain-containing protein 1 (OGFOD1)

Product Code	CSB-EP847635HU-B
Abbreviation	OGFOD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N543
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNGKRPAEPG PARVGKKGKK EVMAEFSDAV TEETLKKQVA EAWSRTPFS HEVIVMDMDP FLHCVIPNFI QSQDFLEGLQ KELMNLDFHE KYNDLYKFQQ SDDLKKRREP HISTLRKILF EDFRSWLSDI SKIDLESTID MSCAKYEFTD ALLCHDDELE GRRIAFILYL VPPWDRSMGG TLDLYSIDEH FQPKQIVKSL IPSWNKLVFF EVSPVSFHQV SEVLSEEKSR LSISGWFHGP SLTRPPNYFE PPIPRSPHIP QDHEILYDWI NPTYLDMDYQ VQIQEEFEES SEILLKEFLK PEKFTKVCEA LEHGHVEWSS RGPPNKRFEY KAEESKLPEI LKECMKLFRS EALFLLLSNF TGLKLHFLAP SEEDEMNDKK EAETTDITEE GTSHSPPEPE NNQMAISNNS QQSNEQTDPE PEENETKKES SVPMCQGELR HWKTGHYTLI HDHСКАEFAL DLILYCGCEG WEPEYGGFTS YIAKGEDEEL LTVNPESNSL ALVYRDRETL KFKVHINHRS LEQKKTFFNR TGFWDFFSFIY YE
Source	E.coli
Target Names	OGFOD1
Protein Names	Recommended name: 2-oxoglutarate and iron-dependent oxygenase domain-containing protein 1 EC= 1.14.11.- Alternative name(s): Termination and polyadenylation 1 homolog
Expression Region	1-542
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.