



Recombinant Human NF-kappa-B-activating protein (NKAP)

Product Code	CSB-MP847637HU
Abbreviation	NKAP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N5F7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAPVSGSRSP DREASGSGGR RRSSSKSPKP SKSARSPRGR RSRSHSCSRS GDRNGLTHQL GGLSQGSRNQ SYRSRERSRS RERPSAPRGI PFASASSSVY YGSYSRPYGS DKPWPSLLDK EREESLRQKR LSERERIGEL GAPEVWGLSP KNPEPDSDEH TPVEDEEPPK TTTSASTSEE EKKKKSSRSK ERSKKRRKKK SSKRKHKKYS ESDSDSDSE TDSSDEDNKR RAKKAKKKEK KKKHRSKKYK KKRSKKSKE SSDSSSKESQ EEFLNPWKD RTKAEPSDL IGPEAPKTLT SQDDKPLNYG HALLPGEAA MAEYVKAGKR IPRRGEIGLT SEEIASFECG GYVMSGSRHR RMEAVRLRKE NQIYSADEKR ALASFNQEER RKRENKILAS FREMVYRKTG GKDDK
Source	Mammalian cell
Target Names	NKAP
Protein Names	Recommended name: NF-kappa-B-activating protein
Expression Region	1-415
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.