



# Recombinant Human Methyltransferase-like protein 13 (METTL13)

<b>Product Code</b>	CSB-EP847648HU
<b>Abbreviation</b>	METTL13
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N6R0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNLLPKSSRE FGSVDYWEKF FQQRGKKAFE WYGTYLELCG VLHKYIKPRE KVLVIGCGNS ELSEQLYDVG YRDIVNIDIS EVVIKQMKEC NATRRPQMSF LKMDMTQMEF PDASFQVVLD KGTLDAVLTD EEEKTLQQVD RMLAEVGRVL QVGGRYLCIS LAQAHILKKA VGHFSREGWM VRVHQQVANSQ DQVLEAEPQF SLPVFAFIMT KFRPVPGSAL QIFELCAQEQ RKPVRLESAE RLAEAVQERQ QYAWLCSQLR RKARLGSVSL DLCDGDTGEP RYTLHVVDSP TVKPSRDNHF AIFIIPQGRE TEWLFMDEG RKQLAASAGF RRLITVALHR GQQYESMDHI QAELSARVME LAPAGMPTQQ QVPFLSVGGD IGVRTVQHGD CSPLSGDYVI EDVQGDDKRY FRRLIFLSNR NVVQSEARLL KDVSHKAQKK RKKDRKKQRP ADAEDLPAAP GQSIDKSYLC CEHHKAMIAG LALLRNPELL LEIPLALLVV GLGGGSLPLF VHDHFPKSCI DAVEIDPSML EVATQWFGFS QSDRMKVHIA DGLDYIASLA GGGEARPCYD VIMFDVDSKD PTLGMSCPPP AFVEQSFLQK VKSILTPEGV FILNLVCRDL GLKDSVLAGL KAVFPLLYVR RIEGEVNEIL FCQLHPEQKL ATPELLETAQ ALERTLRKPG RGWDDTYVLS DMLKTVKIV
<b>Source</b>	E.coli
<b>Target Names</b>	METTL13
<b>Protein Names</b>	Recommended name: Methyltransferase-like protein 13 EC= 2.1.1.-
<b>Expression Region</b>	1-699
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.