



# Recombinant Human RNA polymerase II transcription factor SIII subunit A3 (TCEB3C)

<b>Product Code</b>	CSB-BP847713HU
<b>Abbreviation</b>	TCEB3C
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NG57
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAAGSTTLRA VGKLQVRLAT KTEPKKLEKY LQKLSALPMT ADILAETGIR KTVKRLRKHQ HVGDFARDLA ARWKKLVLDV RNTGPDQPDP EESASRQRFQ EALQEREKAW GFPENATAPR SPSHSPEHRR TARRTPPGQQ RPHPRSPSRE PRAERKRPRM APADSGPHRD PPTRTAPLPM PEGPEPAVPG EQPGRGHAHA AQGGPLLGGQ CQGQPQGEAV GSHSKGHKSS RGASAQKSPV VQESQSERLQ AAGADSAGPK TVPSHVFSEL WDPSEAWMQA NYDLLSAFEA MTSQANPEAL SAPALQEEAA FPGRRVNAKM PVYSGSRPAC QLQVPTLRQQ CLRVPRNNPD ALGDVEGVPI SVLEPVLEGW TPDQLYRTEK DNAALARETD ELWRIHCLQD FKEEKPEHE SWRELYLRLR DAREQRLRVV TTKIRSAREN KPSGRQTKMI CFNSVAKTPY DASRRQEKSA GAADPGNGEM EPAPKPAQSS QAPSGLDGDG GGSVSGGGSS NRHAAPADKT RKQAACKVAP LMAKAIRDYK GRFSRR
<b>Source</b>	Baculovirus
<b>Target Names</b>	ELOA3
<b>Protein Names</b>	Recommended name: RNA polymerase II transcription factor SIII subunit A3 Alternative name(s): Elongin-A3 Short name= EloA3 Transcription elongation factor B polypeptide 3C
<b>Expression Region</b>	1-546
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.