



Recombinant Human Uncharacterized protein C2orf63 (C2orf63)

Product Code	CSB-EP847749HU
Abbreviation	C2orf63
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NHS4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSVHQIRKHA VLPPIICRSD KEFLESVQRY IITETERLGC SEEGPADEYY IYRNVFDKV IEHITAYKSI LTSIKKEYDA FIETIKKDRR TTFCLHGKLK GLAAEPTALV YYRKRTIQLE AKMRIIESNS SKIQSQIDHI KQCRAEYDTK EVKYCTFSKD PSKPIPGMTL QESMNLALT KYMKHLEDKY AEIKQAMLIK YVPAQRKADL DEEMIVLLKR RDVAENLNKK LQFCHQRLQI ISQALSSWVK SDMSSPFQDF VEQIQKTKYL QGDQGIVEEL MEDDPRAKE AEIMLHYIER FNELISLGEY EKAACYAANS PRRILRNIGT MNTFKAVGKI RGKPLPLLLF FEALFITSHA FPCPVDAALT LEGIKCGLSE KRLDLVTNWV TQERLTFSEE AGDVICDYGE QDTYNKAKCL ALAQIVYSEC GLHKKAILCL CKQGQTHRVM EYIQQLKDFD TDDLQLLMS CPQVELIQCL TKELNEKQPS LSFGLAILHL FSADMKKVGI KLLQEINKGG IDAVESLMIN DSFCSIEKWQ EVANICSQNG FDKLSNDITS ILRSQAAVTE ISEEDDAVNL MEHVFW
Source	E.coli
Target Names	CLHC1
Protein Names	Recommended name: Uncharacterized protein C2orf63
Expression Region	1-586
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.