



# Recombinant Human Uncharacterized protein C2orf63 (C2orf63)

<b>Product Code</b>	CSB-MP847749HU
<b>Abbreviation</b>	C2orf63
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NHS4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSVHQIRKHA VLPPIICRSD KEFLESVQRY IITETERLGC SEEGPADEYY IYRNVFDKV IEHITAYKSI LTSIKKEYDA FIETIKKDRR TTFCLHGKLK GLAAEPTALV YYRKRTIQLE AKMRIIESNS SKIQSQIDHI KQCRAEYDTK EVKYCTFSKD PSKPIPGMTL QESMNLALT KYMKHLEDKY AEIKQAMLIK YVPAQRKADL DEEMIVLLKR RDVAENLNKK LQFCHQRLQI ISQALSSWVK SDMSSPFQDF VEQIQKTKYL QGDQGIVEEL MEDDPRAKE AEIMLHYIER FNELISLGEY EKAACYAANS PRRILRNIGT MNTFKAVGKI RGKPLPLLLF FEALFITSHA FPCPVDAALT LEGIKCGLSE KRLDLVTNWV TQERLTFSEE AGDVICDYGE QDTYNKAKCL ALAQIVYSEC GLHKKAILCL CKQGQTHRVM EYIQQKDFD TDDLLQLLMS CPQVELIQCL TKELNEKQPS LSFGLAILHL FSADMKKVGI KLLQEINKGG IDAVESLMIN DSFCSIEKWQ EVANICSQNG FDKLSNDITS ILRSQAAVTE ISEEDDAVNL MEHVFW
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CLHC1
<b>Protein Names</b>	Recommended name: Uncharacterized protein C2orf63
<b>Expression Region</b>	1-586
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.