



# Recombinant *Cricetulus griseus* Guanine nucleotide-binding protein G (s) subunit alpha (GNAS)

<b>Product Code</b>	CSB-MP848127DXU
<b>Abbreviation</b>	GNAS
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8R4A8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Cricetulus griseus</i> (Chinese hamster) ( <i>Cricetulus barabensis griseus</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGCLGNSKTE DQRNEEKAQR EANKKIEKQL QKDKQVYRAT HRLLLLGGAGE SGKSTIVKQM RILHVNGFNG EGGEEDPQAA RSNSDGEKAT KVQDIKNNLK EAIETIVAAM SNLVPPVELA NPENQFRVDY ILSVMNVPNF DFPPEFYEHA KALWEDEGVR ACYERSNEYQ LIDCAQYFLD KIDVVKQADY VPSDQDLLRC RVLTSIGIFET KFQVDKVNFB MFDVGGQRDE RRRKWIQCFND VTAIFVVAS SSYNNMVIRED NQTNRLQEAL NLFKSIWNNR WLRTISVILF LNKQDLLAEK VLAGKSKIED YFPEFARYTT PEDATPEPGE DPRVTRAKYF IRDEFRLIST ASGDGRHYCY PHFTCAVDTE NIRRVFND CR DIIQRMHLRQ YELL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	GNAS
<b>Protein Names</b>	Recommended name: Guanine nucleotide-binding protein G(s) subunit alpha Alternative name(s): Adenylate cyclase-stimulating G alpha protein
<b>Expression Region</b>	1-394
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.