



Recombinant Encephalitozoon cuniculi Probable cell division protein kinase ECU11_1290 (ECU11_1290)

Product Code	CSB-EP848320EKH-B
Abbreviation	ECU11_1290
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8SQU8
Product Type	Recombinant Protein
Immunogen Species	Encephalitozoon cuniculi (strain GB-M1) (Microsporidian parasite)
Purity	>85% (SDS-PAGE)
Sequence	MRS GGRCNSI ENYEKVCRI SSGSFGNVYRV RRKTDNRVFA LKRMNPSMCY DTNGFSILYI REVMILKHIR HRNIMEIEEV VEGCEINDFF IVMECCDTDL RSVIHSVGI GMKAARFLTC QMLKGLKFLH GAGIVHRDLK PSNILLMRDG GLRIADFGLA RAIESQMTNL VVTLWYRPIE ILLGSETYDE SIDMWSVGCV VGEMLRGEPI LAGEGEMDQL DRIFRLLGYP TDADFGLDL PHFKNIRRPS TFEASFEGDF ECYGEEAASF VRNLLSFDPR KRCTASQGLC SGFVADAEEC PGELVDIVGR CTGDI
Source	E.coli
Target Names	ECU11_1290
Protein Names	Recommended name: Probable cell division protein kinase ECU11_1290 EC= 2.7.11.22
Expression Region	1-315
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.