



Recombinant Drosophila melanogaster Protein singed wings 2 (swi2)

Product Code	CSB-YP848367DLU
Abbreviation	swi2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8SXT3
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	A DDSHETSTLL ATITTITTEI IPMPDSGRSC FGWSSAPEAA GGNCRSRNRN GTLKCYGGMN NLAALNQSGK LSRVQPALEM LLCGWPKDGL NHFRDLQKLP RLRSLTIEYS GFTEFKDFDP EMLELHTINI SWTNLSYISS RTFKRVHPLK VLDLRWNQLI QLDGPLLLPR NFEQLYLAGN PWNCTRNFKW LLLQPEKGRL VVDRDELICT DRKYKERQML MVMHYKLELK RQCQSHEDLR NCTCLMHHIL PKTHIPLYTV NCSHLQFHRL PDFLPDNTTT LVINDNMISD INPLRDNPY RHVDMQLEN NQISNVDNLE DTYWLQNFRL LNLRGNNLRK LHVYALDNAL DDNENANLLL LSRNPWHCTC KFGSRMRELL TKYKDIVRDA WNVSCTYRLD DDQLLAKVLT LSRQEMCNLS VEGGTQIHPI DWLNGVLASL IFLILGKLAY DYKYKYYGR VPWIVMKMP
Source	Yeast
Target Names	swi2
Protein Names	Recommended name: Protein singed wings 2
Expression Region	30-499
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.