



Recombinant Drosophila melanogaster Probable E3 ubiquitin-protein ligase sinah (sinah)

Product Code	CSB-EP848376DLU-B
Abbreviation	sinah
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8T3Y0
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSVRNSRPQL SWPERVSPQR TIDTPTASGE MLTRRQSAPA LVVPPEETTH VVVVKRQSPD AAAAGELVPS RRKDSVAVQS GIVATGPLDT TRSGARDDFL MALLECPVCF GYIMPPIMQC PRGHLCSTC RSKLTICPVC RVFMTNIRSL AMEKVASKLI FPCKHSHFGC RARLSYAEKT KHEEDCECRP YFCPYPDDKC SWQGPLRDVY QHLMSSHENV ITMEGNDIIF LATNVNLEGA LDWTMVQSCH GRHFLLSLEK INLGEDCQQY FTACRMIGSM KDAAEFVYNI SLEAYNRTLRL WQSKPRSIRES NFSSFTNADF LVLNKHTVEL FSEDGNLALN VVIRKVEERT N
Source	E.coli
Target Names	sinah
Protein Names	Recommended name: Probable E3 ubiquitin-protein ligase sinah EC= 6.3.2.- Alternative name(s): Sina homolog
Expression Region	1-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.