



Recombinant Human Anterior gradient protein 3 homolog (AGR3)

Product Code	CSB-BP848404HU
Abbreviation	AGR3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TD06
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	IAIKKEKRP PQTLSRGWGD DITWVQTYEE GLFYAQKSKK PLMVIHHLED CQYSQALKKV FAQNEEQEM AQNKFIMLNL MHETTDKNLS PDGQYVPRIM FVDPSLTVRA DIAGRYSNRL YTYEPRDLPL LIENMKKALR LIQSEL
Source	Baculovirus
Target Names	AGR3
Protein Names	Recommended name: Anterior gradient protein 3 homolog Short name= AG-3 Short name= AG3 Short name= hAG-3 Alternative name(s): Breast cancer membrane protein 11
Expression Region	22-166
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.