



Recombinant Human TRAF3-interacting protein 1 (TRAF3IP1)

Product Code	CSB-BP848411HU
Abbreviation	TRAF3IP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TDR0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNAAVVRRRTQ EALGKVIRRP PLTEKLLSKP PFRYLHDIIT EVIRMTGFMK GLYTD AEMKS DNVKDKDAKI SFLQKAIDVV VMVSGEPLLA KPARIVAGHE PERTNELLQI IGCCLNKLKLS SDDAVRRVLA GEKGEVKGRA SLTSRSQELD NKNVREEESR VHKNTEDRGD AEIKERSTSR DRKQKEELKE DRKPREKDKD KEKAKENGGN RHREGERERA KARARPDNER QKDRGNRERD RDSERKKETE RKSEGGKEKE RLRDRDRERD RDKGKDRDRR RVKNGEHSWD LDREKNREHD KPEKKSASSG EMSKKLSDGT FKDSKAETET EISTRASKSL TTKTSKRRSK NSVEGRKEDN ISAKSLDSIV SGINNEPNQE TTTSEIGTKE ANINSTSISD DNSASLRCEN IQPNPTEKQK GDSTSDAEGD AGPAGQDKSE VPETPEIPNE LSSNIRRI PR PGSARPAPPR VKRQDSMEAL QMDRSGSGKT VSNVITESHN SDNEEDDQFV VEAAPQLSEM SEIEMVTAVE LEEEKHGGL VKKILETKKD YEKLQQSPKP GEKERSLFES AWKKEKDIVS KEIEKLRTSI QTLCKSALPL GKIMDYIQED VDAMQNELQM WHSENQRHAE ALQQEQRITD CAVEPLKAEL AELEQLIKDQ QDKICAVKAN ILKNEEKIQK MVYSINLTSR R
Source	Baculovirus
Target Names	TRAF3IP1
Protein Names	Recommended name: TRAF3-interacting protein 1 Alternative name(s): Interleukin-13 receptor alpha 1-binding protein 1 Microtubule-interacting protein associated with TRAF3 Short name= MIP-T3
Expression Region	1-691
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.