



# Recombinant Human TRAF3-interacting protein 1 (TRAF3IP1)

<b>Product Code</b>	CSB-EP848411HU
<b>Abbreviation</b>	TRAF3IP1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8TDR0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MNAAVVRRRTQ EALGKVIRRP PLTEKLLSKP PFRYLHDIIT EVIRMTGFMK GLYTDAEMKS DNVKDKDAKI SFLQKAIDVV VMVSGEPLLA KPARIVAGHE PERTNELLQI IGTKCLNKLS SDDAVRRVLA GEKGEVKGRA SLTSRSQELD NKNVREEESR VHKNTEDRGD AEIKERSTSR DRKQKEELKE DRKPREKDKD KEKAKENGGN RHREGERERA KARARPDNER QKDRGNRERD RDSERKKETE RKSEGGKEKE RLRDRDRERD RDKGKDRDRR RVKNGEHSWD LDREKNREHD KPEKKSASSG EMSKKLSDGT FKDSKAETET EISTRASKSL TTKTSKRRSK NSVEGRKEDN ISAKSLDSIV SGINNEPNQE TTTSEIGTKE ANINSTSISD DNSASLRCEN IQPNPTEKQK GDSTSDAEGD AGPAGQDKSE VPETPEIPNE LSSNIRRI PR PGSARPAPPR VKRQDSMEAL QMDRSGSGKT VSNVITESHN SDNEEDDQFV VEAAPQLSEM SEIEMVTAVE LEEEKHGGL VKKILETKKD YEKLQQSPKP GEKERSLFES AWKKEKDIVS KEIEKLRTSI QTLCKSALPL GKIMDYIQED VDAMQNELQM WHSNRQHAE ALQQEQRITD CAVEPLKAEL AELEQLIKDQ QDKICAVKAN ILKNEEKIQK MVYSINLTSR R
<b>Source</b>	E.coli
<b>Target Names</b>	TRAF3IP1
<b>Protein Names</b>	Recommended name: TRAF3-interacting protein 1 Alternative name(s): Interleukin-13 receptor alpha 1-binding protein 1 Microtubule-interacting protein associated with TRAF3 Short name= MIP-T3
<b>Expression Region</b>	1-691
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.