



Recombinant Human Acid phosphatase-like protein 2 (ACPL2)

Product Code	CSB-EP848416HU
Abbreviation	ACPL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TE99
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	LQFFHLI PVSTPKNGMS SKSRKRIMPD PVTEPPVTDP VYEALLYCNI PSVAERSMEG HAPHHFKLVS VHVFIHGDY YPLYVIPKTK RPEIDCTLVA NRKPYHPKLE AFISHMSKGS GASFESPLNS LPLYPNHPLC EMGELTQTGV VQHLQNGQLL RDIYLKHKHL LPNDWSADQL YLETTGKSRT LQSGLLALYG FLPDFDWKKI YFRHQPSALF CSGSCYCPVR NQYLEKEQRR QYLLRLKNSQ LEKTYGEMAK IVDVPTKQLR AANPIDSMMLC HFCHNVSFPC TRNGCVDMEH FKVIKTHQIE DERERREKKL YFGYSLLGAH PILNQTIGRM QRATEGRKEE LFALYSAHDV TLSPVLSALG LSEARFPRFA ARLIFELWQD REKPSEHSVR ILYNGVDVTF HTSFCQDHHK RSPKPMCPLE NLVRFVKRDM FVALGGSGTN YYDACHREGF
Source	E.coli
Target Names	PXYLP1
Protein Names	Recommended name: Acid phosphatase-like protein 2 EC= 3.1.3.2
Expression Region	24-480
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.