



Recombinant Viral replication protein A28-like (A30L), partial, Yeast

| | |
|--------------------------|---|
| Product Code | CSB-YP848659MHV |
| Abbreviation | A30L |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8V4U9 |
| Product Type | Recombinant Protein |
| Immunogen Species | Monkeypox virus (strain Zaire-96-I-16) (MPX) |
| Purity | >85% (SDS-PAGE) |
| Source | Yeast |
| Target Names | A30L |
| Protein Names | Recommended name: Viral replication protein A28-like Alternative name(s): Protein A30 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | N-terminal His-tagged Tag-Free |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |