



Recombinant Mouse Nitrilase homolog 1 (Nit1)

Product Code	CSB-EP848685MO-B
Abbreviation	Nit1
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8VDK1
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MLGFITRPPH QLLCTGYRLL RTPVLCTQPR PRTMSSSTSW ELPLVAVCQV TSTPNKQENF KTCAELVQEA ARLGACLAFL PEAFDFIARN PAETLLLSEP LNGDLLGQYS QLARECGIWL SLGGFHERGQ DWEQNQKIYN CHVLLNSKGS VVASYRKTHL CDVEIPGQGP MRESNYTKPG GTLEPPVKTP AGKVGLAICY DMRFPELSLK LAQAGAEILT YPSAFGSVTG PAHWEVLLRA RAIESQCYVI AAAQCGRHHE TRASYGHSMV VDPWGTVVAR CSEGPGCLCA RIDLHFLQQM RQHLPVFQHR RPDLYGSLGH PLS
Source	E.coli
Target Names	Nit1
Protein Names	Recommended name: Nitrilase homolog 1 EC= 3.5.-.-
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>