



Recombinant Mouse Protein naked cuticle homolog 2 (Nkd2)

Product Code	CSB-BP848690MO
Abbreviation	Nkd2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VE28
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	GKFQSKHAA AACKRRESPE GDSFVVPAYG SGRRGAEETD RRAGSGVEHR SRDKQELLNG DPKEGPFWDD KGSLEVVLP EKSEGHEGQG QLFSTDDGEK AASREGPLRL SKKHLNIDAL QCDVSVEEDN RQEWFTFLYD FDNSGKVTRE DMSSLMHTIY EVVDASVNHS SGSSKTLRVK LTVSPEPSSK KECPLTGQDR EPTRGRTEIE LTDEPRVADR RLSAYSRKPN ADPQPCSVRV PYCVDENTER RNHYLDLAGI ENYTSKFGPG SPPEQARQEH HGRATHIPSR SRSQESDAHA IHRRSQVLA EHVIPANEP TRALAAQPRI KGQEKQFLRS PKGPGKPLGT PGSGKPGKAL SYCLQAVPLP QSAQDGHHP QPPPQPPQP YGHKRYRQKA REGHSPLKGH GQPTMVEHEV VRDLPPMLGP EGYVMPVVQR HEHHHHHEHH HHHHHHHFHP S
Source	Baculovirus
Target Names	Nkd2
Protein Names	Recommended name: Protein naked cuticle homolog 2 Short name= Naked-2 Short name= mNkd2
Expression Region	2-461
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.