



Recombinant Human Cytosolic arginine sensor for mTORC1 subunit 1 (CASTOR1)

Product Code	CSB-EP848815HU
Abbreviation	GATSL3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WTX7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MELHILEHRV RVLVSVARPG L WLYTHPLIKL LFLPRRSRCK FFSLTETPED YTLMVDEEGF KELPPSEFLQ VAEATWLVLN VSSHSGAAVQ AAGVTKIARS VIAPLAHHV SVLMLSTYQT DFILVREQDL SVVIHTLAQE FDIYREVGGE PVPVTRDDSS NGFPRTQHGP SPTVHPIQSP QNRFCVLTLD PETLPAIATT LIDVLFYSHS TPKEAASSSP EPSSITFFAF SLIEGYISIV MDAETQKKFP SDLLLTSSSG ELWRMVRIGG QPLGFDECGI VAQIAGPLAA ADISAYYIST FNFDHALVPE DGIGSVIEVL QRRQEGLAS
Source	E.coli
Target Names	CASTOR1
Protein Names	Recommended name: GATS-like protein 3
Expression Region	1-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.