



Recombinant Human TBC1 domain family member 22A (TBC1D22A)

Product Code	CSB-MP848818HU
Abbreviation	TBC1D22A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WUA7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MASDGARKQF WKRSNSKLPG SIQHVGGAQH PPFDPPLLHGT LLRSTAKMPT TPVKAKRVST FQEFESNTSD AWDAGEDDDE LLAMAAESLN SEVVMETANR VLRNHSQRQG RPTLQEGPGL QQKPRPEAEP PSPPSGDLRL VKSVSESHTS CPAESASDAA PLQRSQSLPH SATVTLGGTS DPSTLSSSAL SEREASRLDK FKQLLAGPNT DLEELRRLSW SGIPKVRPM TWKLLSGYLP ANVDRRPATL QRKQKEYFAF IEHYYDSRND EVHQDTYRQI HIDIPRMSPE ALILQPKVTE IFERILFIWA IRHPASGYVQ GINDLVTPFF VVFICEYIEA EEVDTVDSVG VPAEVL CNIE ADTYWCMSKL LDGIQDNYTF AQPQIQMKVK MLEELVSRID EQVHRHLDQH EVRYLQFAFR WMNLLMREV PLRCTIRLWD TYQSEPDGFS HFHLYVCAAF LVRWRKEILE EKDFQELLLF LQNLPTAHWD DEDISLLAE AYRLKFAFAD APNHYYK
Source	Mammalian cell
Target Names	TBC1D22A
Protein Names	Recommended name: TBC1 domain family member 22A
Expression Region	1-517
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.