



Recombinant Human Oncoprotein-induced transcript 3 protein (OIT3)

Product Code	CSB-MP848834HU
Abbreviation	OIT3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WWZ8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	DPCSAYIS LNEPWRNTDH QLDESQGPPL CDNHVNGEWY HFTGMAGDAM PTFCIPENHC GTHAPVWLNG SHPLEGDGIV QRQACASFNG NCCLWNTTVE VKACPGGYV YRLTKPSVCF HVYCGHFYDI CDEDCHGSCS DTSECTCAPG TVLGPDRQTC FDENECEQNN GGCSEICVNL KNSYRCECGV GRVLRSDGKT CEDVEGCHNN NGGCSHSCLG SEKGYQCECP RGLVLSEDNH TCQVPVLCKS NAIEVNIPRE LVGGLELFLT NTSCRGVSNV THVNILFSLK TCGTVVDVVN DKIVASNLVT GLPKQTPGSS GDFIIRTSKL LIPVTCEFPR LYTISEGYVP NLRNSPLEIM SRNHGIFPFT LEIFKDNEFE EPYREALPTL KLRDSLYFGI EPVVHVSGLE SLVESCFAFP TSKIDEVLKY YLIRDGCVSD DSVKQYTSRD HLAKHVFQVPV FKFVVGKDHKE VFLHCRVLVC GVLDERSRCA QGCHRRMRRG AGGEDSAGLQ GQTLTGPIR IDWED
Source	Mammalian cell
Target Names	OIT3
Protein Names	Recommended name: Oncoprotein-induced transcript 3 protein Alternative name(s): Liver-specific zona pellucida domain-containing protein
Expression Region	23-545
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.