



# Recombinant Human E3 ubiquitin-protein ligase rififylin (RFFL)

<b>Product Code</b>	CSB-EP848846HU-B
<b>Abbreviation</b>	RFFL
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8WZ73
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MWATCCNWFC LDGQPPEVPP PQGARMQAYS NPGYSSFPSP TGLEPCKSC GAHFANTARK QTCLDCKKNF CMTCSSQVGN GPRCLLLCQR FRATAFQREE LMKMKVKDLR DYLSLHDIST EMCREKEELV LLVLGQQPVI SQEDRTRAST LSPDFPEQQA FLTQPHSSMV PPTSPNLPSS SAQATSVPPA QVQENQQANG HVSQDQEEPV YLESVARVPA EDETQSIDSE DSFVPGRRAS LSDLTDLEDI EGLTVRQLKE ILARNFVNYK GCCEKWELME RVTRLYKDQK GLQHLVSGAE DQNGGAVPSG LEENLCKICM DSPIDCVLLE CGHMTCTKC GKRMNECPIC RQYVIRAVHV FRS
<b>Source</b>	E.coli
<b>Target Names</b>	RFFL
<b>Protein Names</b>	Recommended name: E3 ubiquitin-protein ligase rififylin EC= 6.3.2.- Alternative name(s): Caspase regulator CARP2 Caspases-8 and -10-associated RING finger protein 2 Short name= CARP-2 FYVE-RING finger protein Sakura Short
<b>Expression Region</b>	1-363
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.