



Recombinant Human E3 ubiquitin-protein ligase rififylin (RFFL)

Product Code	CSB-MP848846HU
Abbreviation	RFFL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WZ73
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MWATCCNWFC LDGQPEEVPP PQGARMQAYS NPGYSSFPSP TGLEPSCSKSC GAHFANTARK QTCLDCKKNF CMTCSSQVGN GPRCLLLCQR FRATAFQREE LMKMKVKDLR DYLSLHDIST EMCREKEELV LLVLGQQPVI SQEDRTRAST LSPDFPEQQA FLTQPHSSMV PPTSPNLPSS SAQATSVPPA QVQENQQANG HVSQDQEEPV YLESVARVPA EDETQSIDSE DSFVPGRRAS LSDLTDLEDI EGLTVRQLKE ILARNFVNYK GCCEKWELME RVTRLYKDQK GLQHLVSGAE DQNGGAVPSG LEENLCKICM DSPIDCVLLE CGHMVTCTKC GKRMNECPIC RQYVIRAVHV FRS
Source	Mammalian cell
Target Names	RFFL
Protein Names	Recommended name: E3 ubiquitin-protein ligase rififylin EC= 6.3.2.- Alternative name(s): Caspase regulator CARP2 Caspases-8 and -10-associated RING finger protein 2 Short name= CARP-2 FYVE-RING finger protein Sakura Short
Expression Region	1-363
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.