



Recombinant Acetylnornithine/succinyldiaminopimelate aminotransferase (argD)

Product Code	CSB-MP848863EOD
Abbreviation	argD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8X4S6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	AIEQTAITR ATFDEVILPI YAPAEFIPVK GQGSRIWDQQ GKEYVDFAGG IAVTALGHCH PALVNALKTQ GETLWHISNV FTNEPALRLG RKLIEATFAE RVVFMNSGTE ANETAFKLAR HYACVRHSPF KTKIIAFHNA FHGRSLFTVS VGGQPKYSDG FGPKPSDIIH VPFNDLHAVK AVMDDHTECAV VVEPIQGEGG VTAATPEFLQ GLRELCDQHQ ALLVFDEVQC GMGRTGDLFA YMHYGVTPDI LTSAKALGGG FPISAMLT TA EIASAFHPGS HGSTYGGNPL ACAVAGAAFD IINTPEVLEG IQAKRQRFVD HLQKIDQQYD VFSDIRGMGL LIGAE LKPQY KGQARDFLYA GAEAGVMVLN AGPDVMRFAP SLVVEDADID EGMQRFAHAV AKVVGA
Source	Mammalian cell
Target Names	argD
Protein Names	Recommended name: Acetylnornithine/succinyldiaminopimelate aminotransferase Short name= ACOAT Short name= DapATase Short name= Succinyldiaminopimelate transferase EC= 2.6.1.11 EC= 2.6.1.17
Expression Region	2-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.