



Recombinant 23S rRNA (uracil (747)-C (5))-methyltransferase RlmC (rlmC)

Product Code	CSB-MP848876EOD
Abbreviation	rlmC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8X6Q5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MQCALYDAGR CRSCQWITQP IPEQLSAKTA DLKNLLADFP VEEWCAPVSG PEQGFRNKAK MVVSGSVEKP LLGMLHRDGT PEDLCDCPLY PASFAPVFAA LKPFARAGL TPYNVARKRG ELKYILLTES QSDGGMMLRF VLRSETKLAQ LRKALPWLQA QLPQLKVITV NIQPVHMAIM EGETEIYLTE QQALVERFND VPLWIRPQSF FQTNPAVASQ LYATARDWVR QLPVKHMWDL FCGVGGFGLH CATPDMQLTG IEIALEAIAC AKQSAAELGL TRLQFQALDS TQFATAQGEV PELVLVNPPR RGIGKPLCDY LSTMAPRFII YSSCNAQTMA KDIRELPGYR IERVQLFDMF PHTAHYEVLTL LLVKM
Source	Mammalian cell
Target Names	rlmC
Protein Names	Recommended name: 23S rRNA (uracil(747)-C(5))-methyltransferase RlmC EC= 2.1.1.189 Alternative name(s): 23S rRNA(m5U747)-methyltransferase
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.