



Recombinant Ribosomal large subunit pseudouridine synthase E (rluE)

Product Code	CSB-BP848879EOD
Abbreviation	rluE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8X724
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	≥85% (SDS-PAGE)
Sequence	MRQFIISNT MQKTSFRNHQ VKRFSSQRST RRPENQPTR VILFNKPYDV LPQFTDEAGR KTLKEFIPVQ GVYAAGRLDR DSEGLLVLTN NGALQARLTQ PGKRTGKIYY VQVEGIPTQD ALEALRNGVT LNDGPTLPAG AELVDEPAWL WPRNPPIRER KSIPTSWLKI TLYEGRNRQV RRMTAHVGFP TLRLIRYAMG DYSLDNLANG EWREATE
Source	Baculovirus
Target Names	rluE
Protein Names	Recommended name: Ribosomal large subunit pseudouridine synthase E EC=5.4.99.20 Alternative name(s): rRNA pseudouridylate synthase E rRNA-uridine isomerase E
Expression Region	1-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.