



Recombinant Putative transposase yncI (yncI)

Product Code	CSB-MP848884EOD
Abbreviation	yncI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8X7R5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MSIQSLLDYI SVIPDIRQQG KVKHKLSDIL FLTVCAVIAG ADEWQEIEDF GHERLEWLKK YGDFDNGIPV DDTIARVVSNI DLSLAFEKIF IEWMQECHEI TDGEIIAIDG KTIRGSFDKG KRKGAIHMVS AFSNENGVVL GQVKTEAKSN EITAIPELLN LLDLKKNLIT IDAMGCQKDI ASKIKDKKAD YLLAVKGNQG KLHHAFEEKF PVNVFSNYKG DSFSTQEISH GRKETRLHIV SNVTPEFCDF EFEWKGLKKL CVALSFRQKK EDKSAEGVSI RYYISSKDMD AKEFAHAIRA HWLIEHSLHW VLDVKMNEDA SRIRRGNAAE IISGIKKMAL NLLRDCKDIK GGVKKRKKV ALNTCYIEEV LASCSELGFR TDKMKNLTQI
Source	Mammalian cell
Target Names	yncI
Protein Names	Recommended name: Putative transposase yncI
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.