



Recombinant Biotin synthase (bioB)

| | |
|--------------------------|---|
| Product Code | CSB-BP848887EOD |
| Abbreviation | bioB |
| Storage | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p> |
| Uniprot No. | Q8X825 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli O157:H7 |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MAHRPRWTL S QVTELF EKPL LDLLFEAQQV HRQHFDPRQV QVSTLLSIKT GACPEDCKYC PQSSRYKTGL EAERLMEVEQ VLESARKAKA AGSTRFCMGA AWKNPHERDM SYLEQMVQGV KAMGLEACMT LGTLESSEQA RLANAGLDYY NHNLDTSPEF YGNIITRTY QERLDTLEKV RDAGIKVCSG GIVGLGETVK DRAGLLLQLA NLPTPPESVP INMLVKVKG T PLADNDDVDA FDFIRTIAVA RIMMPTS YVR LSAGREQMNE QTQAMCFMAG ANSIFYGCKL LTPNPEEDK DLQLFRKLGL NPQQTAVLAG DNEQQQRLEQ ALMTPDTDEY YNAAAL |
| Source | Baculovirus |
| Target Names | bioB |
| Protein Names | Recommended name: Biotin synthase EC= 2.8.1.6 |
| Expression Region | 1-346 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p> |