



Recombinant GMP reductase (guaC)

Product Code	CSB-MP848897EOD
Abbreviation	guaC
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8X984
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	<pre>MRIEDLKLKLG FKDVLIRPKR STLKSRSDVE LERQFTFKHS GQSWSGVPPII AANMDTVGTF SMASALASFD ILTAVHKHYS VEEWQAFINN SSADVLKHVM VSTGTSDADF EKTKQILDNL PALNFVCIDV ANGYSEHFVQ FVAKAREAWP TKTICAGNVV TGEMCEELIL SGADIVKVGI GPGSVCTTRV KTGVGYPQLS AVIECADA AH GLGGMIISDG GCTTPGDVAK AFGGGADFVM LGGMLAGHEE SGGRIVEENG EKFMLFYGMS SESAMKRHVG GVAEYRAAEG KTVKLPLRGP VENTARDILG GLRSACTYVG ASRLKELTKR TTFIRVQEQE NRIFNNL</pre>
Source	Mammalian cell
Target Names	guaC
Protein Names	Recommended name: GMP reductase EC= 1.7.1.7 Alternative name(s): Guanosine 5'-monophosphate oxidoreductase Short name= Guanosine monophosphate reductase
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>