



Recombinant UPF0214 protein yfeW (yfeW)

Product Code	CSB-EP848918EOD-B
Abbreviation	yfeW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8XBJ0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	A KYPVLTESSP EKAGFNVERL NQMDRWISQQ VDAGYPGVNL LIIKDNQIVY RKAWGAAKKY DGSVLMAQPV KATTGTLYDL ASNTKMYATN FALQKLMSEG KLHPDDLIK YISGFADSPN DTIKGKNTLR ISDLLHHS GG FPADPQYPNK AVAGALYSQD KGQTLEMIKR TPLEYQPGSK HIYSDVDYML LGFIVESVTG QPLDRYVEES IYRPLGLTHT VFNPLLKGFK PQQIAATELN GNTRDGVVHF PNIRTSTHWG QVHDEKAFYS MGGVSGHAGL FSNTGDI AVL MQTMLNGGGY GDVQLFSAET VKMFTTSSKE DATFGLGWRV NGNATMTPTF GTLASPQTYG HTGWTGTVTV IDPVNHMAIV MLSNKPSPV ADPQKNPNMF ESGQLPIATY GWVVDQVYAA LKQK
Source	E.coli
Target Names	yfeW
Protein Names	Recommended name: UPF0214 protein yfeW
Expression Region	20-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.