



Recombinant Lysine--tRNA ligase (lysS)

Product Code	CSB-MP848933EOD
Abbreviation	lysS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8XD57
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	SEQHAQGAD AVVDLNNELK TRREKLANLR EQGIAFPNDF RRDHTSDQLH AEFDGKENEE LEALNIEVAV AGRMMTRRIM GKASFVTLQD VGGRIQLYVA RDDLPEGVYN EQFKKWDLGD ILGAKGKLFK TKTGELSIHC TELRLLTKAL RPLPDKFHGL QDQEARYRQR YLDLISNDES RNTFKVRSQI LSGIRQFMVN RGFMEVETPM MQVIPGGAAA RPFITHNAL DLDMYLRIAP ELYLKRLVVG GFERVFEINR NFRNEGISVR HNPEFTMMEL YMAYADYKDL IELTESLFRT LAQDILGKTE VTYGDVTLDF GKPFEKLTMR EAIKKYRPET DMADLDNFDS AKAIVESIGI HVEKSWGLGR IVTEIFEEVA EAHLIQPTFI TEYPAEVSPL ARRNDVNPEI TDRFEFFIGG REIGNGFSEL NDAEDQAQRF LDQVAAKDAG DDEAMFYDED YVTALEHGLP PTAGLGIGID RMVMLFTNSH TIRDVILFPA MRPVK
Source	Mammalian cell
Target Names	lysS
Protein Names	Recommended name: Lysine--tRNA ligase EC= 6.1.1.6 Alternative name(s): Lysyl-tRNA synthetase Short name= LysRS
Expression Region	2-505
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.