



Recombinant Human Homeobox protein Hox-D11 (HOXD11)

Product Code	CSB-MP010683HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P31277
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNDFDECGQS AASMYLPGCA YYVAPSDFAS KPSFLSQPSS CQMTFPYSSN LAPHVQPVRE VAFRDYGLER AKWPYRGGGG GGSAGGGSSG GGPGGGGGGA GGYAPYYAAA AAAAAAAAAA EEAAMQRELL PPAGRRPDVL FKAPEPVCAA PGPPHGPAGA ASNFYSAVGR NGILPQGFDQ FYEAAPGPPF AGPQPPPPPA PPQPEGAADK GDPRTGAGGG GGSPCTKATP GSEPKGAAEG SGGDGEGPPG EAGA EKSSSA VAPQRSRKKR CPYTKYQIRE LEREFFFNVY INKEKRLQLS RMLNLTDRQV KIWFQNRMMK EKKLNRDRLQ YFTGNPLF
Source	Mammalian cell
Target Names	HOXD11
Protein Names	Recommended name: Homeobox protein Hox-D11 Alternative name(s): Homeobox protein Hox-4F
Expression Region	1-338
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene belongs to the homeobox family of genes. The homeobox genes encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, located on different chromosomes, consisting of 9 to 11 genes arranged in tandem. This gene is one of several homeobox HOXD genes located in a cluster on chromosome 2. Deletions that remove the entire HOXD gene cluster or the 5 end of this cluster have been associated with severe limb and genital abnormalities. The product of the mouse Hoxd11 gene plays a role in forelimb morphogenesis.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.