



Recombinant *Pyrobaculum aerophilum* 50S ribosomal protein L22P (rpl22p)

Product Code	CSB-BP849497FHU
Abbreviation	rpl22p
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8ZWH7
Product Type	Recombinant Protein
Immunogen Species	<i>Pyrobaculum aerophilum</i> (strain ATCC 51768 / IM2 / DSM 7523 / JCM 9630 / NBRC 100827)
Purity	>85% (SDS-PAGE)
Sequence	MPRHQSYSLS DEDIIELVFR NYGVKITPEQ IARAYAPEQR MSWKKSVEVA RFIKGMTLKQ AKTWLEDVVK LKRPVPIKTF KKKQAHHATP WEGWPVAKWP VKVARRFLQV LENLENNARF RGLDVDRVVI VHAAAHKGFK IPNIMPRAFG RATRFDEQTV NVELVAVELP KEVLPKRYKL NLVKR
Source	Baculovirus
Target Names	rpl22
Protein Names	Recommended name: 50S ribosomal protein L22P
Expression Region	1-185
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.