



# Recombinant Chicken Heart- and neural crest derivatives-expressed protein 1 (HAND1)

<b>Product Code</b>	CSB-MP849534CH
<b>Abbreviation</b>	HAND1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q90691
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MNLVGGYQHH HHHHHHHHML HEPFLFGPAA RCHQERAYFP GWVLNPAEVT PELHGQSPAY GPAEFGSGGA GRLEALSGRL GRRKGVGGPK KERRRTESIN SAFAELRECI PNVPADTKLS KIKTLRLATS YIAYLMEVLA RDSQPGEPEG FKAELKKADG RENKRKRETQ PEVYSQPLAH GEKCLKGRTG WPQQVWALEL NP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	HAND1
<b>Protein Names</b>	Recommended name: Heart- and neural crest derivatives-expressed protein 1 Alternative name(s): Extraembryonic tissues, heart, autonomic nervous system and neural crest derivatives-expressed protein 1 Short name= eHAND
<b>Expression Region</b>	1-202
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein belongs to the basic helix-loop-helix family of transcription factors. This gene product is one of two closely related family members, the HAND proteins, which are asymmetrically expressed in the developing ventricular chambers and play an essential role in cardiac morphogenesis. Working in a complementary fashion, they function in the formation of the right ventricle and aortic arch arteries, implicating them as mediators of congenital heart disease. In addition, it has been suggested that this transcription factor may be required for early trophoblast differentiation.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.