



Recombinant Maize streak virus genotype D Replication-associated protein (C1/C2)

Product Code	CSB-MP849636MIAF
Abbreviation	C1/C2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91MF8
Product Type	Recombinant Protein
Immunogen Species	Maize streak virus genotype D (isolate Raw) (MSV)
Purity	>85% (SDS-PAGE)
Sequence	MTSSSSNRQF LHRTANTFLT YPQCPEHPEI ISQRIWDLVG RWNPLYIICA QEAHEDGNMH LHALIQTDKQ VRTTDSRFFD IDGFHPNIQS AMSPNKVRDY ILKEPLALFE RGTFVPRKKT FLGNSSKGNS EKKPSKDEIM QDIISHATSK PEYLSMVRKS FPYDWATKLQ YFEYSANKLF PDIQEEFINP HPTSEPDLLC NESIKDWLQP NIYQQADNGS RKQSLYIVGP TRTGKSTWAR SLGRHNYWQN NVDWSSYNED AIYNIVDDIP FKYCPCWKQL VGCQKEFVVN PKYGKKKKVQ MKSKPTIILA NSDEDWMKEM TPGQLEYFEA NCMIIYMSPG EKWYSPPQLP PTEEV
Source	Mammalian cell
Target Names	C1/C2
Protein Names	Recommended name: Replication-associated protein Short name= Rep EC= 2.7.7.- EC= 3.1.21.-
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.