



Recombinant Mouse Cold shock domain-containing protein E1 (Csde1)

Product Code	CSB-EP849654MO-B
Abbreviation	Csde1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91W50
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MSFDPNLLHN NGHNGYPNGT SAALRETGVI EKLLTSYGFI QCSEQRARLF FHCSQYNGNL QDLKVGDDVE FEVSSDRRTG KPIAIKLVKI KPEIHPEERM NGQVVCAPPH NLESKSPAAP GQSPTGSVCY ERNGEVFYLT YTSEDVEGNV QLETGDKINF VIDNNKHTGA VSARNIMLLK KKQARCQGVV CAMKEAFGFI ERGDVVKEIF FHYSEFKGDL ETLQPGDDVE FTIKDRNGKE VATDVRLLPQ GTVIFEDISI EHFEGTVTKV IPKVPSKNQN DPLPGRIKVD FVIPKELPFG DKDTKSKVTL LEGDHVRFNI STDRRDKLER ATNIEVLSNT FQFTNEAREM GVIAAMRDGF GFIKCVDRDA RMFFHFSEIL DGNQLHIAD EFTTVVPDML SAQRNHAIRI KKLPGKTVSF HSHSDHRFLG TVEKEATFSN PKTTSPNKGK DKEAEDGIIA YDDCGVKLTI AFQAKDVEGS TSPQIGDKVE FSISDKQRPG QQIATCVRLR GRNSNSKRLL GYVATLKDNF GFJETANHDK EIFFHYSEFS GDVDSLELGD MVEYSLSKGG GNVKSAEKNV KAHSVNGITE EANTPIYSGK VIRPLRGVDP TQIEYQGMIE IVEEGDMKGE VYPFGVIGMA NKGDCQLQGE SVKFQLCVLG QNAQTMAYNI TPLRRATVEC VKDQFGFINY EVGDSKKLFF HVKEVQDGVE LQAGDEVEFS VILNQRGTGKC SACNVWRVCE GPKAVAAPRP DRLVNRLKNI TLDDASAPRL MVLQRQPRGPD NSMGFGAERK IRQAGVID
Source	E.coli
Target Names	Csde1
Protein Names	Recommended name: Cold shock domain-containing protein E1
Expression Region	1-798
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.