



Recombinant Arabidopsis thaliana Probable WRKY transcription factor 69 (WRKY69)

Product Code	CSB-BP850034DOA
Abbreviation	WRKY69
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q93WV5
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MHRRAAIQES DDEEDETYND VVPESPSSCE DSKISKPTPK KRRNVEKRVV SVPIADVEGS KSRGEVYPPS DSWAWRKYGQ KPIKGSYPR GYYRCSSSKG CPARKQVERS RVDPSKLMIT YACDHNHPFP SSSANTKSHH RSSVVLKTAK KEEYEEEEEE ELTVTAAEEP PAGLDLSHVD SPLLLGGCYS EIGFEGWFYD ASISSSSGSS NFLDVTLEGR FSVGQEEDES LFGDLGLPD CASVFRRGTV ATEEQHRRCD FGAIPFCDSS R
Source	Baculovirus
Target Names	WRKY69
Protein Names	Recommended name: Probable WRKY transcription factor 69 Alternative name(s): WRKY DNA-binding protein 69
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.