



Recombinant Human cAMP-dependent protein kinase type I-beta regulatory subunit (PRKAR1B)

Product Code	CSB-EP018696HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P31321
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ASPPACPSE EDESLKGCEL YVQLHGIQQV LKDCIVHLCI SKPERPMKFL REHFEEKLEKE ENRQILARQK SNSQSDSHDE EVSPTPPNPV VKARRRRGGV SAEVYTEEDA VSYVRKVIPK DYKTMTALAK AISKNVLFAH LDDNERSDIF DAMFPVTHIA GETVIQQGNE GDNFYVVDQG EVDVYVNGEW VTNISEGGSF GELALIYGTP RAATVKAKTD LKLWGIDRDS YRRILMGSTL RKRKMYEEFL SKVSILESLE KWERLTVADA LEPVQFEDGE KIVVQGE PGD DFYIITEGTA SVLQRRSPNE EYVEVGR LGP SDYFGEIALL LNRPRAATVV ARGPLKCVKL DRPRFERVLG PCSEILKRNI QRYNSFISLT V
Source	E.coli
Target Names	PRKAR1B
Protein Names	Recommended name: cAMP-dependent protein kinase type I-beta regulatory subunit
Expression Region	2-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.